



NEW 200 EVO

cGMP compact Part Washer

iwt[®]
a **TECNIPLAST** company
PHARMA

NEW!

200 EVO

cGMP compact Part Washer

Engineered to clean product contact parts from solid, semi-solid, liquid and aseptic manufacturing sites, from tableting machines to sterile filling lines

GAMP5, cGMP, ASME-BPE, ISPE baseline, FDA (including 21CFR part 11), CE Marking and UL Listing compliant

- ▶ **Modular design** to allow quick turnaround projects
- ▶ **Compact solution** to simplify installations
- ▶ **Multiple configurations available**
- ▶ **Three different sizes:** 600, 1000, 1400 liters (160, 265, 370 gal)
- ▶ **Up to 30 fully customizable washing programs** from pre-washing to drying and cooling
- ▶ **PanelPC configuration for remote connectivity and INDUSTRY 4.0 integration;** cutting edge automation system GAMP5
- ▶ **HEPA filtered drying system**
- ▶ **Multiple water inlet and chemical dosing systems**
- ▶ **Single Pass Rinse feature**
- ▶ **Full in-house design and manufacturing of customized racks** connected via automatic quick-lock patented system
- ▶ **Comprehensive set of handling solutions** for manoeuvring, stocking and loading/unloading of custom racks, inclusive of height-adjustable trolleys to enhance operator's ergonomics

200 EVO MODEL	060	100	140
Chamber Nominal Volume [m ³ – ft ³]	> 0.6 – 21.2	> 1.00 – 35.3	> 1.4 – 49.4
Footprint [m ² – ft ²]	~ 1.2 – 12.9	~ 1.8 – 19.4	~ 2.5 – 26.9
Unit height [mm - ft]	2500 - 8.2		
Unit total volume [m ³ – ft ³]	3 – 105.8	4.5 – 159.1	6.3 – 220.6
Ratio	20%	22%	22%

Best-in-class ratio between chamber and total volumes in the smallest footprint for parts cleaning



to find out more information:

info@iwtpharma.com | www.iwtpharma.com

CAT MKTBCR000130L- Rev Nov 2022

IWT reserves the right to change or modify product and/or specifications without notice or obligation.