



CLASSIC line

Article Number: 29-200667-030

Product Range

Safety storage cabinet

Product Category

Type 90 - flammable liquids

General

Line	CLASSIC
Version	standard
Size	M
Type	90
Storage type	passive storage
Functionality	Storage
Earthing connection	Yes
Connection socket	Yes

Certification

TÜV Süd certificate no.	Z1A 012906 0520
TÜV-Süd certificate type	GS
TÜV-Süd High Quality Zertifikat	Z2 012906 0513
DIN EN 14470-1	Yes
DIN EN 16121/16122	Yes

Technical specifications

Outer width	594 mm
Outer depth	612 mm
Outer height	2045 mm
Max. load bearing capacity per storage area	75 kg
Max. load	360 kg
Max. total weight	635 kg
Max. volume of single container	10 L
Collection volume	11 L
Weight (product)	265 kg

Cabinet

Material of cabinet cladding	Sheet steel
Cabinet colour	RAL 7035
Material of inner lining	Decor panels
Colour of inner lining	RAL 7035
Base	Classical
Base height	89 mm
Base colour	RAL 7015

Doors

No. of doors	1 Stk
Kind of door	Wing door
Material of door cover	Sheet steel
Door colour	RAL 1018
Door hinge	Right
Door operation	Classical
Door technology	Classical
Door handle	Bow handle matt nickel
Lock	Cylinder lock
Locking system	Open
Locking points	1-point

Interior fittings

Material of interior fittings	Sheet steel
Colour of interior fittings	RAL 7035
Bottom tray	Fixed
Storage shelf	Yes
Storage levels	4 Stk
Storage areas	4 Stk
setting height	394 mm